

# CS4

FURNACE APPARATUS  
FOR KLY5-A AND MFK1-FA  
KAPPABRIDGES

# AGICO

ADVANCED GEOSCIENCE INSTRUMENTS COMPANY

Supplement for KLY5-A and MFK1-FA Kappabridges for measurements of high-temperature variations of magnetic susceptibility.

## General Description

The **CS4 High Temperature Furnace Apparatus** is optional attachment for **KLY5-A** and **MFK1-FA/A** Kappabridges designed for measurement of the temperature variations of low-field magnetic susceptibility of minerals, rocks and synthetic materials in the temperature range from **ambient temperature to 700°C**. Measurements can be performed under the protective argon atmosphere to prevent oxidation of measured specimen.

The apparatus consists of a non-magnetic furnace with a special platinum temperature sensor, a temperature control unit, and a cooling water reservoir. The specimen (up to 0.25 cm<sup>3</sup> in volume) is placed in a silica glass vessel, heated by a platinum wire, and the temperature is measured by the temperature sensor.

The quasicontinuous measurement process is fully automated, being controlled by the software **Su-fyfe**. Data processing software **Cureval** serves for advanced analysis of thermomagnetic curves such as empty furnace measurement subtraction, Curie temperature estimation and separation of ferromagnetic and paramagnetic part of susceptibility.

## CS4 Furnace Apparatus Comprising

- CS4 Temperature Control Unit
- Furnace
- Water Cooling Reservoir
- Temperature Sensor
- Argon Gas Flowmeter
- Power Supply Unit
- Specimen Vessels
- Set of Interconnecting Cables
- Measuring Software, CUREVAL Software
- User's Manual

## Main Features

- Lowest detectable susceptibility change  $1 \times 10^{-7}$  SI.
- Measurement at high temperatures up to 700°C.
- Software controlled heating and cooling modes.
- Measurement in the air or argon atmosphere.
- 500 to 700 pairs of susceptibility and temperature measurements to define a thermomagnetic curve.

## Technical specifications

Nominal specimen volume:	0.25 cm <sup>3</sup>
Inner diameter of measuring vessel:	6 mm
Temperature range:	room temp up to 700°C
Accuracy of temperature determination:	$\pm 2^\circ\text{C}$
Sensitivity to susceptibility changes:	$1 \times 10^{-7}$ SI
Power requirements:	100 - 240 V, 50/60 Hz, 350 VA



LABORATORY INSTRUMENTS FOR MEASUREMENT OF MAGNETIC PROPERTIES OF ROCKS

# CS-L

CRYOSTAT APPARATUS  
FOR KLY SERIES AND MFK1-FA / A  
KAPPABRIDGES

# AGICO

ADVANCED GEOSCIENCE INSTRUMENTS COMPANY

Supplement for KLY5-A and MFK1-FA Kappabridges for measurements of low-temperature variations of magnetic susceptibility.

## General Description

The **CS-L Low Temperature Cryostat Apparatus** is optional attachment for **KLY series** and **MFK1-FA/A** Kappabridges designed for measurement of the temperature variation of low-field magnetic susceptibility of minerals, rocks and synthetic materials in the temperature range from **-192°C to ambient temperature**. CS4 (or CS-3) apparatus is mandatory prerequisite for CS-L system.

The specimen is placed in a measuring vessel which is cooled inside the cryostat by liquid nitrogen and then heated spontaneously to a given temperature. The argon gas is needed for depleting the liquid nitrogen out of cryostat.

The quasi-continuous measurement process is fully automated, being controlled by the software **Su-fyte**. Data processing software **Cureval** serves for advanced analysis of thermomagnetic curves such as empty furnace measurement subtraction, Curie temperature estimation and separation of ferromagnetic and paramagnetic part of susceptibility.

## CS-L Cryostat Apparatus Comprising

Cryostat

Temperature Sensor

Specimen Vessels

Pot and Funnel for Liquid Nitrogen

Argon Blow System

Measuring Software

## Main Features

Measurement at low temperatures from -192° C

Cooled by liquid nitrogen

Estimation of ratio between ferromagnetic and paramagnetic part of susceptibility

## Technical specifications

Nominal specimen volume: 0.25 cm<sup>3</sup>

Inner diameter of measuring vessel: 6 mm

Temperature range: -192°C - ambient temperature

Accuracy of temperature determination: ±2 °C

Sensitivity to susceptibility changes: 1 x 10<sup>-7</sup> SI

Power requirements: 100 - 240 V, 50/60 Hz, 350 VA



## AGICO, Inc.

Advanced Geoscience Instruments Company  
Ječná 29a, CZ – 621 00 Brno, Czech Republic

Tel.: +420 511 116 303 Fax.: +420 541 634 328 E-mail: [agico@agico.cz](mailto:agico@agico.cz) Web: [www.agico.com](http://www.agico.com)